

# Calculation of Biochemical Methane Potential (BMP)\*

Sasha D. Hafner  
sasha.hafner@eng.au.dk

March 13, 2020

*Document number 200. File version 1.3. This document is from  
the Standard BMP Methods collection.*<sup>†</sup>

## 1 Background

This document describes how to calculate biochemical methane potential (also called biomethane potential) (BMP) from laboratory measurements made in a batch BMP test. The calculations are based on standardized CH<sub>4</sub> volume produced in bottles with: 1) inoculum only and 2) with substrate and inoculum, along with information on the quantity of inoculum and substrate volatile solids (VS) added to each bottle. For details on calculating CH<sub>4</sub> production, see method-specific documents from the Standard BMP Methods collection<sup>1</sup>.

### 1.1 Selection of a BMP duration

The time at which to evaluate BMP, i.e., the length of the incubation, should be at least as long as the 1% net duration<sup>2</sup>. Regardless, it is important that the time is identical for both inoculum-only and inoculum + substrate bottles when carrying out calculation of BMP. This does not mean it cannot vary among substrates from within the same BMP test.

---

\*Recommended citation: Hafner, S.D. Calculation of Biochemical Methane Potential (BMP). Standard BMP Methods document 201, version 1.4. Available online: <https://www.dbfz.de/en/BMP> (accessed on Feb 28, 2020).

Or see <https://www.dbfz.de/en/BMP> for a BibTeX file that can be imported into citation management software.

<sup>†</sup>For more information and other documents, visit <https://www.dbfz.de/en/BMP>. For document version history or to propose changes, visit <https://github.com/sashahafner/BMP-methods>.

<sup>1</sup> Available methods are: volumetric (document 201), manometric (document 202), gravimetric (document 203), and gas density (document 204), all of which can be downloaded from <https://www.dbfz.de/en/BMP>.

<sup>2</sup> For a detailed description, see document 100 from <https://www.dbfz.de/en/BMP>.

## 1.2 SMP and BMP

The term specific methane production (SMP) is used to refer to  $\text{CH}_4$  yield from a particular substrate, in the same units as BMP<sup>3</sup>. BMP is in fact simply SMP at a single suitable duration. SMP curves, or SMP for all measurement intervals, show the development of SMP over time, reflecting kinetics of a particular test, and are commonly included in BMP reports. The calculations described below are used for both SMP and BMP.

## 2 Calculation of SMP and BMP

These calculations require the following variables. Units may differ, but typical units are listed below.

- $V_{\text{CH}_4, S, i, t}$ , the standardized volume of  $\text{CH}_4$  produced in bottle  $i$  containing inoculum and substrate at time  $t$  (mL)
- $V_{\text{CH}_4, I, j, t}$ , the standardized volume of  $\text{CH}_4$  produced in bottle  $j$  with inoculum only (“blank” bottle) at time  $t$
- $m_{I, i}$ , the mass of inoculum (typically as-measured (fresh) mass) originally added to bottle  $i$
- $m_{VS, S, i}$ , the mass of substrate volatile solids (VS) originally added to bottle  $i$
- $n$ , the number of replicate bottles with inoculum + substrate (typically 3, the minimum)
- $k$ , the number of replicate bottles with inoculum only (“blanks”, typically 3, the minimum)

Methane productivity of inoculum is calculated separately for each inoculum-only bottle at each measurement time using Eq. (1).

$$v_{\text{CH}_4, I, j, t} = V_{\text{CH}_4, I, j, t} / m_{I, j} \quad (1)$$

And from these, a mean value is calculated as

$$\bar{v}_{\text{CH}_4, I, t} = \sum_{j=1}^k v_{\text{CH}_4, I, j, t} / k \quad (2)$$

where  $k$  = the number of inoculum-only bottles.

---

<sup>3</sup> standardized  $\text{CH}_4$  volume (dry, 0°C, 101.325 kPa, referred to as “normal” volume) per unit mass of substrate VS added (often written as  $\text{NmL}_{\text{CH}_4} \text{ gVS}^{-1}$ ).

Net CH<sub>4</sub> production from inoculum + substrate bottles, i.e., an estimate of CH<sub>4</sub> production derived from substrate only<sup>4</sup> is calculated as given in Eq. (3).

$$V_{CH_4,S,i,net} = V_{CH_4,S,i} - \bar{v}_{CH_4,I} \cdot m_{I,i} \quad (3)$$

Note that the units on inoculum mass are completely irrelevant and have no effect of results, as long as they are sufficiently precise. Fresh (wet) mass is recommended, although dry or VS mass could be used.<sup>5</sup> Equation (3) may be simplified

Bottle yield is calculated by normalizing net CH<sub>4</sub> production by substrate VS mass:

$$B_i = V_{CH_4,S,i,net} / m_{VS,S,i} \quad (4)$$

Finally, BMP for a particular substrate is taken as the mean of these values.

$$\bar{B} = \sum_{i=1}^n B_i / n \quad (5)$$

where  $n$  is the number of replicate bottles.

### 3 Calculation of random error

Calculation of random error in BMP estimates must include two significant sources: variation in apparent yield among substrate bottles, and variation in apparent inoculum yield in inoculum-only bottles. Uncertainty in determination of substrate VS may sometimes be a significant source of error, and, while not required, its inclusion is recommended. Here, these all are quantified using standard error and will be referred to as  $s_{\bar{x},1}$  (substrate yield),  $s_{\bar{x},2}$  (inoculum yield), and  $s_{\bar{x},3}$  (substrate VS). Note that  $s_{\bar{x},1}$  and  $s_{\bar{x},2}$  may include many sources of error that collectively contribute to the observed value. Other sources of random error are assumed to be small: determination of inoculum and substrate mass in particular. Systematic error, which may be more important in some cases, is not included here.

Units on all three standard errors are the units of the final BMP estimates, e.g., standardized CH<sub>4</sub> volume from substrate in mL per g substrate VS. They can be added together to provide a total estimate with:

$$s_{\bar{x},BMP} = \sqrt{\sum_{i=1}^3 s_{\bar{x},i}^2} \quad (6)$$

---

<sup>4</sup>This calculation is based on the assumption of additivity for CH<sub>4</sub> production, i.e., production of CH<sub>4</sub> from inoculum is not affected by the presence of substrate. This is almost certainly incorrect, but similar results even when varying the inoculum-to-substrate ratio suggest it is not a large source of error.

<sup>5</sup>Any error in determination of inoculum dry matter or VS here is exactly canceled by the combination of Eqs. (1) and (3), so has no effect.

Given  $s_{\bar{x},BMP}$ , standard deviation  $s_{BMP}$  can be calculated by multiplying by  $\sqrt{n}$ . Although this approach seems to be preferred in the literature, the interpretation is a bit ambiguous.

The value of  $s_{\bar{x},1}$  is calculated from yield values calculated for individual bottles:

$$s_{\bar{x},1} = \sqrt{\sum_{i=1}^n \frac{(B_i - \bar{B})^2}{n-1}} \quad (7)$$

To calculate  $s_{\bar{x},1}$  the standard error of normalized inoculum-only  $\text{CH}_4$  production (e.g., volume of  $\text{CH}_4$  in mL per g inoculum mass) is first calculated from individual values determined from each individual inoculum-only bottle.

$$s_{\bar{x},I} = \sqrt{\sum_{j=1}^k \frac{(v_{CH_4,I,j} - \bar{v}_{CH_4,I})^2}{k-1}} \quad (8)$$

For each individual substrate bottle an estimate of  $s_{\bar{x},1}$  is made as:

$$s_{\bar{x},2,i} = s_{\bar{x},I} \cdot m_{I,i} \quad (9)$$

Equation (9) includes an assumption that error in inoculum mass determination is negligible compared to error in inoculum  $\text{CH}_4$  yield, which is reasonable. Given these values,  $s_{\bar{x},2}$  is then calculated with:

$$s_{\bar{x},2} = \sqrt{\frac{\sum_{i=1}^n s_{\bar{x},2,i}^2}{n}} \quad (10)$$

The third source of random error is determined from  $s_{\bar{x},VS,i}$ , which is the standard error in substrate VS mass addition, calculated as the product of substrate (fresh) mass (assumed to include only negligible random error) and substrate VS concentration standard error, which is assumed to be the only significant source of random error for  $s_{\bar{x},3}$ . Resulting units for  $s_{\bar{x},VS,i}$  are the same as the substrate VS mass units. Given this value, an estimate of  $s_{\bar{x},3}$  can be made for each bottle with:

$$s_{\bar{x},3,i} = \frac{s_{\bar{x},VS,i}}{m_{VS,i}} \cdot B_i \quad (11)$$

Equation (11) differs from Eq. (9) because for  $s_{\bar{x},3}$  error is propagated through division and not subtraction. (Note that standard error of substrate VS concentration could be used directly in Eq. (11).) Given these values for individual bottles,  $s_{\bar{x},3}$  is calculated from Eq. (12).

$$s_{\bar{x},3} = \sqrt{\frac{\sum_{i=1}^n s_{\bar{x},3,i}^2}{n}} \quad (12)$$